

ESR 22-21698 Study Requisition Form

Patient ID: _____

| Collection Day | Test | Tube Type | Collection Time (24 Hour Clock) | Initials & Date Collected |
|--------------------------------|---|---|------------------------------------|------------------------------|
| Screening | COVID-19 | COVID-19 Test | ____:____ | |
| Day 1 | SARS-CoV-2 RBD-IgG Antibody Level | (1) 6.0 ml Gold Top SST | ____:____ | |
| Day 30 – No Friday Collections | T-Cell Assay & AZD7442 Serum Concentration | (3) 8.5 ml Yellow Top ACD, (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |
| Day 60 | AZD7442 Serum Concentration | (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |
| Day 90 | SARS-CoV-2 RBD-IgG Antibody Level & AZD7442 Serum Concentration | (1) 6.0 ml Gold Top SST, (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |
| Day 120 | AZD7442 Serum Concentration | (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |
| Day 150 | AZD7442 Serum Concentration | (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |
| Day 180 | SARS-CoV-2 RBD-IgG Antibody Level & AZD7442 Serum Concentration | (1) 6.0 ml Gold Top SST, (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |
| Day 210 | AZD7442 Serum Concentration | (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |
| Day 270 | AZD7442 Serum Concentration | (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |
| Day 360 | AZD7442 Serum Concentration | (1) 3.5 ml Gold Top SST, (2) Cryovials | ____:____ | |